

 ATHENAGTX**PS1201**

PS1201 PRESSURE SLEEVE INFUSION DEVICE

The *PS1201* is a small, light weight, autonomous, electronic alternative to traditional manual pressure IV bags. It automatically applies and maintains constant pressure during IV infusion. A “set-it and forget-it” device that frees providers’ hands. *PS1201* is ideal for interosseous infusions requiring sustained higher pressures than conventional IV treatments. Highly portable and quickly deployed, *PS1201* streamlines both patient and provider movement.

Performance under pressure.



Size: *With sleeve wrapped around case - 6" X 2.7" X 1.5"*
(15.25cm X 6.8cm X 3.8cm)

Open - 6" X 11.8" X 1.25"
(15.25cm X 30cm X 3.2cm)

Weight: 6 oz (170 gm) with battery

Control: On/Off rocker switch activates internal circuit with pressure transducer

Indicator: LED - Green = operating, Orange = battery low

Power: One CR123 lithium battery

Battery life: 3 hours, approx.

Sleeve capacity: One 500 ml or One Liter IV bag, Integral hanging hook for bag

Nominal Pressure: 300mmHg, manual push button vent.

 sales@athenagtx.com

 515.288.3360

 athenagtx.com

 ATHENAGTX

5900 NW 86TH STREET, SUITE 300, JOHNSTON, IA 50131

650-0036-01 Rev A