



Quick and easy deployment



Rugged construction for any field application

TacVac™ PORTABLE SUCTION DEVICE

TacVac offers all the capability and performance of traditional battery powered suction devices in a more mobile package. Initially designed for far-forward military field applications, TacVac features light-weight, rugged construction and a very small form factor. TacVac can be collapsed into an ultra-compact storage configuration and instantly deployed when needed.



FDA cleared

Operational configuration



The TacVac is infinitely carry-able and instantly deployable.

- Quick assembly
- Strong suction power
- Low Cost

Size: 7.25" x 3.0" x 3.5"

Weight: 9.2 ounces w/
two CR123 A batteries

Batteries: Two CR123A
Lithium/Manganese Dioxide
(Li/MnO₂) - commercially
available

**Collection Container
Capacity:** 500mL

Accessories: fluid collection
container and suction tip
(quantities sold separately)

*Suction motor has the
smallest form factor of
anything in it's class*



*Compact, lightweight
and portable*



APPLICATIONS

- EMS
- TACTICAL
- SEARCH AND RESCUE
- SKILLED NURSING
- HOME HEALTHCARE

*Ergonomic
design for easy
handling*



CONTACT US TODAY TO GET MORE INFORMATION OR TO SCHEDULE A DEMO.



sales@athenagtx.com



515.288.3360



athenagtx.com