

INFECTIOUS DISEASE TRANSPORT



iDT INFECTIOUS DISEASE TRANSPORT

MONITOR BOTH PATIENTS AND PROVIDERS



Wirelessly monitor and transmit vital signs data for multiple patients and care providers through isolation and PPE

Isolation and personal protective equipment places special requirements on BOTH patients and providers. *The Wireless Vital Signs Monitor (WVSM®) wirelessly transmits patient data through isolation equipment.* This allows clinicians to monitor both the patient and care providers on a hand held device or smart phone from an area outside the isolation zone. *SurgeCase* can charge up to 10 WVSMs and provide monitoring for up to 10 patients, including remote monitoring of multiple patients simultaneously.

The **WVSM** is small, light and easily worn by a patient in isolation as well as providers under protective gear. Incident commanders can monitor both patients and providers from a central monitoring device even during transport.




 AUTOMATICALLY MONITORS SPO2, HR, NIBP & ECG
 
 WIRELESS MONITORING RANGE OF 200 YARDS
 
 REMOTE CONNECTIVITY USING ADMS



**INFECTIOUS
DISEASE
TRANSPORT**

The Infectious Disease transport kit with WVSMs includes everything providers need to monitor patients and clinicians during an incident

- *Up to Ten (10) WVSM[®]s with peripherals (Measuring NIBP, HR, SpO2 and ECG)*
- *Optional centralized monitoring device*

WVSM FEATURES

The small and portable *Wireless Vital Signs Monitor (WVSM[®])* collects the data you need to make critical decisions, even in remote, dangerous and multi-patient environments.

- *Automatically monitors SpO2, HR, NIBP and ECG*
- *Weighs just one pound*
- *A fraction of the cost of a larger single-patient monitor*
- *Automatically records and transmits 4.5 hours of patient history*
- *Wireless range of up to 200 yards*
- *Up to 7.5 hours of battery life*
- *Color-coded vitals provide the data medics need without overwhelming first responders*
- *Stays with the patient from point-of-injury until treatment*



Optional ruggedized Surge Case for charging and storage of 10 WVSMs plus accessories



Wireless monitoring of up to 20 patients on your compatible device



Contact us today to get more info or to schedule a demo.



SALES@ATHENAGTX.COM



515.288.3360



ATHENAGTX.COM